## LDL Cholesterol

Analytical Performance Specification

RCV Calculation

Matrix	BV Estimate	median CV estimate	lower CI limit	higher CI limit	Comments	Date Updated
Serum/plasma	•	8.3	6.6	10.1		21/9/2023
Serum/plasma	**	26.6	18.5	37.6		21/9/2023

## Colour Key for References

Included In Meta Analysis Not Included I		s						
Reference	Dataset	Estimate of CVI	Estimate of CVG	Gender	Age	State of Well Being	Matrix	Sampling Interval
Age dependence of within-subject biological variation of nine common clinical chemistry analytes Carobene A, Graziani MS, Lo Cascio C, Tretti L, Cremonese E, Yabarek T, Gambaro G and Ceriotti F., Clin Chem Lab Med, 841-4, 50, 2012	View	11.2		Mixed		Healthy	Serum	1 per Week
View								
Short-term intraindividual variability in lipoprotein measurements: the atherosclerosis risk communities (ARIC) study Chambless LE, McMahon RP, Brown SA, Patsch W, Heiss G and Shen Y, Am J Epidem, 1069-1081, 0136, 1992 View	View	6.4		Mixed		Healthy	Plasma	1 per Week
Short-term intraindividual variability in lipoprotein measurements: the atherosclerosis risk communities (ARIC) study	View	7.8		Mixed		Healthy	Plasma	1 per Week

Brown SA, Patsch W, Heiss G and Shen Y, Am J Epidem,

Chambless LE, McMahon RP,

1069-1081, 0136, 1992

View

Variability of plasma lipids and lipoproteins: the Jerusalem Lipid

and seradarem Eipia			EFLM Bio	ogical Variatio	n			
Research Clinic Study Friedlander Y, Kark JD and Stein Y, Clin Chem, 1121-6, 31, 1985	View	12.8		Males	40.0	Healthy	Plasma	1 per Month
View								
Variability of plasma lipids and lipoproteins: the Jerusalem Lipid Research Clinic Study Friedlander Y, Kark JD and Stein Y, Clin Chem, 1121-6, 31, 1985	View	12.5		Males	60.0	Healthy	Plasma	1 per Month
View								
Significance of various parameters derived from biological variability for lipid and lipoprotein analyses Pagani F and Panteghini M, Clin Biochem, 415-20, 26, 1993	View	5.6	21.3	Mixed		Healthy	Serum	1 per Week
View								
Variability in cholesterol measuraments: comparison of calculated and direct LDL cholesterol determinations  Schectman G, Patsches M and Sasse EA, Clin Chem, 732-7,	View	6.6		Males		Healthy	Serum	1 per Month
42, 1996								
View								
Variación biológica de lípids y lipoproteinas séricas en indivíduos sanos y pacientes con enfermedad renal crónica terminal  Corte Z and Venta R, Rev Lab Clin, 88-95, 7, 2014	View	7.3	43.2	Mixed		Healthy	Serum	1 per Week
View								
Variation in concentrations of lipids, lipoprotein lipids, and apolipoproteins A-I and B in plasma from healthy women	View	13.6	18.5	Females		Healthy	Plasma	1 per
Brown SA, Boerwinkle E, Kashanian FK, Swanson N and Patsch W, Clin Chem, 207-10,						·		Other

30, IYYU			EFLM Biol	logical Variatio	n			
View								
Biological variation of glucose, insulin and lípids in lean, insulinsensitie and obese, insulin-resistant chinese males without diabetis Rizi EP, Baig S, Toh SA and Khoo CM, Pathology, 510-2, 48, 2016	View	9.8	19.8	Males		Healthy	Serum	1 per Week
View								
Age dependence of within-subject biological variation of nine common clinical chemistry analytes Carobene A, Graziani MS, Lo Cascio C, Tretti L, Cremonese E, Yabarek T, Gambaro G and Ceriotti F., Clin Chem Lab Med, 841-4, 50, 2012	View	8.4		Mixed		Healthy	Serum	1 per Week
View								
Biological variation of glucose, insulin and lípids in lean, insulinsensitie and obese, insulin-resistant chinese males without diabetis  Rizi EP, Baig S, Toh SA and Khoo CM, Pathology, 510-2,	View	10.5	16.1	Males		Healthy	Serum	1 per Week
48, 2016								
View								
A longitudinal study of the biological variability of plasma lipoproteins in healthy young adults Mjøs O, Rao S,Bjøru, Henden T, Thelle D,Førde O and Miller N, Atherosclerosis, 75-81, 34, 1979	View	10.4		Males	28.0	Healthy	Serum	1 per Other
View								
The EuBIVAS: withinand between-subject biological variation data for electrolytes, lipids, urea, uric acid, total protein, total bilirubin, direct bilirubin and glucose  Aarsand AK, Diaz-Garzón J,	View	8.46	28.2	Mixed		Healthy	Serum	1 per Week

Fernandez-Calle P, et al, Clin Chem, 1380-93, 64(9), 2018						
View						
Biological variability of Cholesterol, Triglyceride, Low- and High-density lipoprotein Cholesterol, Lipoprotein(a), and Apolipoproteins A-I and B Marcovina SM, Gaur VP and Albers JJ, Clin Chem, 574-8, 40, 1994	View	9.2	Mixed	Healthy	Plasma	2 per Month
View						
Biological variation of lipids and lipoproteins in fingerstick blood Kafonek SD, Donovan L, Lovejoy KL and Bachorik PS, Clin Chem, 2002-7, 42, 1996	View	7.5	Mixed 40.3	Healthy	Serum	1 per Month
View						
Biological variation of lipids and lipoproteins in fingerstick blood Kafonek SD, Donovan L, Lovejoy KL and Bachorik PS, Clin Chem, 2002-7, 42, 1996	View	5.4	Mixed 40.3	Healthy	Whole blood	1 per Month
View						
Within-person fluctuations of serum cholesterol and lipoproteins Mogadam M, Ahmed S, Mensch A and Godwin I, Arch Intern Med, 1645-8, 150, 1990	View	9.6	Mixed	Healthy	Serum	1 per Week
View						
Diurnal and monthly intra-individual variability of the concentration of lipids, lipoproteins and apoproteins Wasenius A, Stugaard M, Otterstad JE, Froyshov D, Scand J Clin Lab Invest, 632-42, 50, 1990	View	5.2	Mixed	Healthy	Serum	1 per Week
View						
A longitudinal study of the biological variability of plasma lipoproteins						

in healthy young adults						Serum	1 per
Mjøs O, Rao S,Bjøru, Henden T, Thelle D,Førde O and Miller N, Atherosclerosis, 75-81, 34, 1979	View	11.8		Females 28.0	Healthy	Serum	Other
View							
Significance of various							
parameters derived							
from biological							
variability for lipid and							1 per
lipoprotein analyses	View	5.1	26.3	Mixed	Healthy	Serum	Week
Pagani F and Panteghini M, Clin Biochem, 415-20, 26, 1993							
View							
Variability in cholesterol							
measuraments:							
comparison of							
calculated and direct							
LDL cholesterol	View	6.5		Males	Healthy	Serum	1 per
determinations							Month
Schectman G, Patsches M and Sasse EA, Clin Chem, 732-7, 42, 1996							
View							
Variación biológica de							
lípids y lipoproteinas							
séricas en indivíduos							
sanos y pacientes con							
enfermedad renal	View	7.0	33.4	Mixed	Healthy	Serum	1 per Week
crónica terminal							vveek
Corte Z and Venta R, Rev Lab Clin, 88-95, 7, 2014							
View							
Biological variation data							
applied to the selection							
of serum lipid ratios							
used as risk markers of						C	1 per
coronary heart disease.	View	13.7	24.5	Males	Healthy	Serum	Month
Ortolá J, Castiñeiras MJ,							
Fuentes-Arderiu X., Clin Chem, 56-9, 38, 1992							
View							
Biological variation of							
glucose, insulin and							
lípids in lean, insulin-							
sensitie and obese,							
insulin-resistant chinese males without diabetis	View	13.9	25.8	Males	Healthy	Serum	1 per Week
Rizi EP, Baig S, Toh SA and  Khoo CM Pathology 510-2							

KITOO CIVI, I attrology , JTO 4, 48, 2016 View Variability of plasma lipids and lipoproteins: the Jerusalem Lipid Research Clinic Study 1 per Plasma 12.1 View Females 40.0 Healthy Month Friedlander Y, Kark JD and Stein Y, Clin Chem, 1121-6, 31, 1985 View Variability of plasma lipids and lipoproteins: the Jerusalem Lipid Research Clinic Study 1 per Plasma 11.3 View Females 60.0 Healthy Month Friedlander Y, Kark JD and Stein Y, Clin Chem, 1121-6, 31, 1985 View Variability of plasma lipids and lipoproteins: the Jerusalem Lipid Research Clinic Study 1 per Plasma 13.7 View Healthy Males 17.0 Month Friedlander Y, Kark JD and Stein Y, Clin Chem, 1121-6, 31, 1985 View The EuBIVAS: withinand between-subject biological variation data for electrolytes, lipids, urea, uric acid, total protein, total bilirubin, 1 per Serum View 6.82 21.6 Females 60.0 Healthy direct bilirubin and Week glucose Aarsand AK, Diaz-Garzón J, Fernández-Calle P, et al, Clin Chem, 1380-93, 64(9), 2018 View Prognostic utility of biochemical markers of cardiovascular risk: impact of biological variability 1 per Serum 10.3 View 19.3 Mixed Healthy Month Alexander KS, Kazmierczak SC, Snyder CK, Oberdorf JA and Farrell DH, Clin Chem Lab Med, 1875-82, 51(9), 2013 View Circadian rhythmic

			EFEINI BIO	ogical Variatior	1			
variations in serum concentrations of clinically important		0.6			20.5		Serum	
lipids Rivera-Coll A, Fuentes-Arderiu X and Díez-Noguera A, Clin Chem, 1549-53, 40, 1994	View	8.6		Mixed	29.5	Healthy		4 per Ho
View								
Prognostic utility of biochemical markers of cardiovascular risk: impact of biological variability Alexander KS, Kazmierczak SC, Snyder CK, Oberdorf JA and Farrell DH, Clin Chem Lab Med , 1875-82, 51(9), 2013	View	7.8	20.4	Mixed		Healthy	Serum	2 per Week
View								
Variación biológica de lípids y lipoproteinas séricas en indivíduos sanos y pacientes con enfermedad renal crónica terminal Corte Z and Venta R, Rev Lab Clin, 88-95, 7, 2014	View	11.9	74.2	Mixed		Non- healthy	Serum	1 per Month
View								
Diurnal and monthly intra-individual variability of the concentration of lipids, lipoproteins and apoproteins Wasenius A, Stugaard M, Otterstad JE, Froyshov D, Scand J Clin Lab Invest, 632-42, 50, 1990	View	5.1		Mixed		Healthy	Serum	3 per Hou
View								
The EuBIVAS: within- and between-subject biological variation data for electrolytes, lipids, urea, uric acid, total protein, total bilirubin, direct bilirubin and glucose Aarsand AK, Diaz-Garzón J, Fernández-Calle P, et al, Clin Chem, 1380-93, 64(9), 2018	View	8.55	27.0	Males	35.0	Healthy	Serum	1 per Week

			EFLM Biolo	ogical Variation				
and hemostatic risk								
indicators in men and								
women.								
Riese H, Vrikkotte TGM, Meijer P, Kluft C, De Geus EJC, Fibrinolysis & Proteolysis, 1- 12, 14, 2001	View	7.6	29.8	Females	30.19	Healthy	Serum	2 per Day
View								
Covariance of metabolic								
and hemostatic risk								
indicators in men and								
women.								
Riese H, Vrikkotte TGM, Meijer P, Kluft C, De Geus EJC, Fibrinolysis & Proteolysis, 1- 12, 14, 2001	View	6.5	21.8	Females	36.71	Healthy	Serum	2 per Day
View								
Covariance of metabolic								
and hemostatic risk								
indicators in men and								
women.								
Riese H, Vrikkotte TGM, Meijer P, Kluft C, De Geus EJC, Fibrinolysis & Proteolysis, 1-	View	6.5	23.5	Males	45.23	Healthy	Serum	2 per Day
12, 14, 2001 View								
Age dependence of								
within-subject								
biological variation of nine common clinical								
chemistry analytes	View	6.6		Mixed		Healthy	Serum	1 per
Carobene A, Graziani MS, Lo Cascio C, Tretti L, Cremonese E, Yabarek T, Gambaro G and Ceriotti F., Clin Chem Lab Med, 841-4, 50, 2012								Week
View								
Biological variation of								
glucose, insulin and								
lípids in lean, insulin-								
sensitie and obese,								
insulin-resistant chinese	View	17.4	32.8	Males		Non-	Serum	1 per
males without diabetis						healthy		Week
Rizi EP, Baig S, Toh SA and Khoo CM, Pathology , 510-2, 48, 2016								
View								
View  Analytical and biological								
Analytical and biological								

## **EFLM Biological Variation** Health Study/Study of Latinos 1 per Thyagarajan B, Howard AG, Serum 27.42 View Mixed Healthy Month Durazo-Arvizu R, Eckfeldt JH, Gellman MD, Kim RS, Liu K, Mendez AJ, Penedo FJ, Talavera GA, Youngblood ME, Zhao L and Sotres-Alvarez D, Clin Chim Acta, 129-137, 463, 2016 View Analytical and biological variability in biomarker measurement in the **Hispanic Community** Health Study/Study of Latinos Thyagarajan B, Howard AG, 1 per Serum 9.7 View Mixed Healthy Month Durazo-Arvizu R, Eckfeldt JH, Gellman MD, Kim RS, Liu K, Mendez AJ, Penedo FJ, Talavera GA, Youngblood ME, Zhao L and Sotres-Alvarez D, Clin Chim Acta, 129-137, 463, 2016 View Biological variability of lipoproteins and apolipoproteins in patients referred to a "No Non-**Plasma** 7.6 lipid clinic View Females 9.0 Sampling healthy Intervals" Kafonek SD, Derby CA and Bachorik PS, Clin Chem, 864-72, 38, 1992 View The EuBIVAS: withinand between-subject biological variation data for electrolytes, lipids, urea, uric acid, total protein, total bilirubin, 1 per Serum 9.24 22.5 View Females 34.0 Healthy Week direct bilirubin and glucose Aarsand AK, Diaz-Garzón J, Fernández-Calle P, et al, Clin Chem, 1380-93, 64(9), 2018 View Biological variability of lipoproteins and apolipoproteins in patients referred to a "No Non-Plasma 8.0 lipid clinic View Females 52.0 Sampling healthy Intervals" Kafonek SD, Derby CA and Bachorik PS, Clin Chem, 864-

View								
Biological variability of lipoproteins and								
apolipoproteins in								
patients referred to a								"No
lipid clinic	View	7.0		Males	51.0	Non-	Plasma	Samplin
Kafonek SD, Derby CA and Bachorik PS, Clin Chem, 864- 72, 38, 1992						healthy		Intervals
View								
Biological variability of								
lipoproteins and								
apolipoproteins in								
patients referred to a						<b>b</b> 1 -		"No
lipid clinic	View	7.7		Males	9.0	Non- healthy	Plasma	Samplin
Kafonek SD, Derby CA and Bachorik PS, Clin Chem, 864- 72, 38, 1992						y		Interval
View								
Biological variability of								
lipoproteins and								
apolipoproteins in								
patients referred to a								"No
lipid clinic	View	7.8		Mixed	28.0	Non- healthy	Plasma	Samplin
Kafonek SD, Derby CA and Bachorik PS, Clin Chem, 864- 72, 38, 1992								Interval
View								
Variability of plasma								
lipids and lipoproteins:								
the Jerusalem Lipid								
Research Clinic Study	View	12.5		Females	17.0	Healthy	Plasma	1 per
Friedlander Y, Kark JD and	V ICVV	14.7		i eiliales	17.0	nearthy		Month
Stein Y, Clin Chem, 1121-6, 31, 1985								
View								
Variación biológica de								
lípids y lipoproteinas								
séricas en indivíduos								
sanos y pacientes con								1
enfermedad renal	View	9.8	53.6	Mixed		Non- healthy	Serum	1 per Month
crónica terminal								
Corte Z and Venta R, Rev Lab Clin, 88-95, 7, 2014								
View								

08/11/2023, 09:33 EFLM Biological Variation



European Federation of Clinical Chemistry and Laboratory Medicine
Via Carlo Farini 81
20159 Milano - Italy

This website and its content is copyright of EFLM. All rights reserved. You may not, except with our express written permission, distribute or commercially exploit the content (contact Silvia Cattaneo, EFLM Office for more information, silvia.cattaneo@eflm.eu). You may, however, print ordownload to a local hard disk extracts for your personal and non-commercial use only. If you present the content you must acknowledge the website as the source of the material. You cannot transmit the content or store it in any other website or other form of electronic retrieval system.

**Disclaimer**